

Heavy-Duty Restraints From Morrison Medical



Heavy-Duty Wrist Restraints

- Comfortable blue neoprene over nylon webbing.
- Velcro® to assist restraint placement.
- Quick release plastic buckles for fast on/off.
- Each side has two straps with double D-ring closures for sturdy two-point anchoring.
- Reusable.
- Latex-Free.
- FDA approved.

1172: Heavy-Duty Wrist Restraint, pair

Please call for pricing.



Heavy-Duty Ankle Restraints

- Comfortable red neoprene over nylon webbing.
- Velcro® to assist restraint placement.
- Quick release plastic buckles for fast on/off.
- Each side has two straps with double D-ring closures for sturdy two-point anchoring.
- Reusable.
- Latex-Free.
- FDA approved.

1173: Heavy-Duty Ankle Restraint, pair

Please call for pricing.

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*