

Morrison Medical's Big Blue Head Immobilizer



Big Blue Head Immobilizer *Please Call for Pricing*

Item #1420:
Big Blue Head Immobilizer, each


Item #1420-6:
*Big Blue Head Immobilizer,
case of 6*

Item #1432:
Pediatric Head Blocks, one pair

- Vinyl coated foam for easy disinfecting.
- Extra large pillows provide maximum support.
- Easy access to ears for patient assessment.
- Pillows adjust and firmly attach to base with wide Velcro®.
- Secures to backboard with mounted straps and D-rings.
- Includes base with attached straps, two large pillows and two head/chin straps.
- Adaptable for pediatric patients with the use of pediatric head blocks (#1432).
- Available in blue, orange and black.

Pediatric Head Immobilizer Blocks

- Compatible with Big Blue and most other popular head immobilizers.
- Vinyl coated foam for easy disinfecting.
- Easy access to ears for patient assessment.
- Pillows adjust and firmly attach to base with Velcro®.
- 2½" x 5½" x 3¾" pillows.
- Available in blue.



*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*