

## Morrison Medical's Head Vise III - With Loop Strap

The Head Vise III, a popular selection from the Morrison Medical disposable head immobilizer line, is now available with an adhesive-free loop head strap. This all weather, latex-free head immobilizer will fit virtually all backboards and patients. The hook engagable base allows the wedges to be repositioned to customize the fit to your patient. The extra wide 3" wedges provide great stability and patient comfort. The adjustable loop strap further stabilizes the patient and secures to the hook portion of the wedges. Each immobilizer comes with: one foam base with pressure sensitive adhesive for application, two 3" wide foam wedges with hook sections for securing to base and securing head strap, and one loop head strap. The Head Vise III is both disposable and affordable. Sold individually or in cases of 24.

### *Head Vise III Please Call for Pricing*

*Item #1471:  
Head Vise III with loop strap, each*

*Item #1471-24:  
Head Vise III with loop strap, case of 24*



- Latex-Free
- Includes one adhesive-free loop head strap.
- Foam base measures 4 x 15".
- Extra wide foam wedges (3") for greater stability.
- Single use, disposable.
- Quick and easy to apply.
- Patient's ears are easily accessible.
- Universal – will fit virtually all backboards and patients.

*For more information, call 800-438-6677,  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).*