


## Disposable Pressure Pad - *Complete with cover to allow for reuse*



- Combination of medical grade foam elastomer (similar to gel)
- Cushions patients against pain and pressure
- Adhesive backing keeps pad in place without mess
- Helps support concavities in patient's spinal area
- In the 99th percentile of humans from upper back to sacrum
- Assists in eliminating patient movement while on board
- Hypoallergenic - no latex as found in most standard foams
- Hydrophobic - designed to NOT stick to wounds
- Assists in pediatric torso to head relationship by ATLS standards
- Molded contour of spinal area improves lower back support
- Safe to use in water rescue
- Measures 14"x36" making it Universal to any backboard
- Single use - may qualify for insurance reimbursement

Item 3041: Disposable Pressure Pad, with cover, each  
*Please Call for Pricing*



**Molded contour for support and stability**



*For more information, call 800-438-6677,  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).*