

## Backboards from Morrison Medical




- Central concavity to center patient and decrease pressure and pain.
- Doubles as a pediatric board with special slotted area that uses standard straps to fit smaller patients. (Pediatric pins available)
- 100% x-ray translucent.
- 450, 600 and 1000 pound load capacity.
- Tapered for easy extrication in confined areas.
- Stores easily in emergency vehicles at 72" x 16" x 1¾".
- Custom monogramming available.
- Available in Red, Orange, Yellow, Blue, Green and Lime Green.
- Competitively priced.

**For more information,  
call 800-438-6677  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).**

*Please call for pricing.*

- #3053 Trauma Board 450, No Pins, Each
- #3054 Trauma Board 450, W/Pins, Each
- #3055 Trauma Board 600, No Pins, Each
- #3056 Trauma Board 600, W/Pins, Each
- #3057 Trauma Board 1000, No Pins, Each
- #3058 Trauma Board 1000, W/Pins, Each

Pedi-pins are available by request.



**Available in 6 vibrant colors!**



## How do Morrison Medical's backboards compare to the competition?

	Morrison Medical Trauma 450	Pro Lite Eco Spineboard	Ferno Najo Lite	LSP SpineX
Weight	13 lbs.	16 lbs.	14.5 lbs.	17 lbs.
Dimensions	72 x 16 x 1.75"	72 x 16 x 2.25"	72 x 16 x 1.75"	72.5 x 18 x 2.5"
Concave	YES	NO	NO	NO
Pediatric built-in	YES	NO	NO	YES
Load capacity	450 lbs.	500 lbs.	450 lbs.	350 lbs.
	Morrison Medical Trauma 600	Ferno Najo Redihold	Iron Duck Ultra Vue	LSP Xtra
Weight	16 lbs.	16 lbs.	16 lbs.	15.25 lbs.
Dimensions	72 x 16 x 1.75"	72 x 16 x 1.75"	72 x 16 x 1.75"	71 x 15.75 x 1.85"
Concave	YES	NO	NO	NO
Pediatric built-in	YES	NO	NO	NO
Load capacity	600 lbs.	600 lbs.	600 lbs.	500 lbs.
	Morrison Medical Trauma 1000	Pro Lite XT Spineboard	Ferno Najo Mega 1000	Bak-Pak Backboard
Weight	13.5 lbs.	14 lbs.	20 lbs.	11.5 lbs.
Dimensions	72 x 16 x 1.75"	72 x 16 x 2.25"	72 x 18 x 2.25"	72 x 16 x 2.25"
Concave	YES	NO	NO	NO
Pediatric built-in	YES	NO	NO	NO
Load capacity	1,000 lbs.	1,000 lbs.	1,000 lbs.	1,000 lbs.