

## Folding Splints

For use in limb immobilization, Morrison Medical folding splints are available in a variety of sizes and foam options. In addition to the options featured below, visit [www.morrisonmed.com](http://www.morrisonmed.com) to see the complete line of splints including leg splints, Cold Splints, Vel Splints, Wood Board Splints, and more.

- Feature 5/8" punch-out holes for tie downs
- Micro-perforated edges to prevent paper cuts
- 4 styles: 3 strips of 1/2" thick foam (each side & center); full foam 3/8"; center strip 3/8" foam; and plain (no foam)
- 24" and 34" lengths offer an extra center crease for firm support of thin limbs
- Available in cases (call 800-438-6677 for details)
- Registered by FDA

### Corrugated Plastic Folding Splints



2507 Full foam, 12"  
2550 Full foam, 18"  
2560 Full foam, 24"  
2575 Full foam, 34"



2548 3 piece foam, 12"  
2558 3 piece foam, 18"  
2568 3 piece foam, 24"  
2578 3 piece foam, 34"



2504 Center foam, 12"  
2530 Center foam, 18"  
2540 Center foam, 24"  
2573 Center foam, 34"



2500 Plain, 12"  
2510 Plain, 18"  
2520 Plain, 24"  
2570 Plain, 34"

### Cardboard Folding Splints



1507 Full foam, 12"  
1550 Full foam, 18"  
1560 Full foam, 24"  
1576 Full foam, 34"



1548 3 piece foam, 12"  
1558 3 piece foam, 18"  
1568 3 piece foam, 24"  
1588 3 piece foam, 34"



1504 Center foam, 12"  
1530 Center foam, 18"  
1540 Center foam, 24"  
1573 Center foam, 34"



1500 Plain, 12"  
1510 Plain, 18"  
1520 Plain, 24"  
1570 Plain, 34"

*For more information, call 800-438-6677,  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).*