

The Hip Hugger™

The Hip Hugger™ System

- Intended for suspected pelvic /hip injury stabilization.
- Eases discomfort in transport.
- Cleans easily.
- Color coded straps for easy identification in hurried situations.

As easy as...



The Hip Huggie™

The Hip Hugger™ System

- Intended for suspected pelvic/hip injury stabilization.
- Eases discomfort in transport.
- Cleans easily.
- Color and number coded straps for easy application.
- Smaller and more economical version of the original Hip Hugger.
- Comes with carrying case.



For optimal fit, tighten straps labeled "1" at the same time, followed by "2" straps.

For more information, call 800-438-6677,
or visit www.MorrisonMed.com.