

## Morrison Medical's Head Vise II



*Head Vise II  
Please Call for Pricing*

*Item #1465:  
Head Vise II, each*

*Item #1465SS:  
Head Vise II with two straps, each*



- Made of sturdy moisture resistant cardboard.
- Includes comfortable foam head pad.
- Quickly attaches to board with EZ peel tabs, no straps or clips needed.
- X-ray translucent.
- Patient's ears are easily accessible.
- Universal and adjustable – will fit virtually all backboards and patients.
- Available with two blue all weather straps, 2 x 36".
- Individually sealed in plastic to be protect from the elements until ready to use.
- Compact size.
- Stores flat.
- Single use, disposable.
- Writable surface.

*For more information, call 800-438-6677,  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).*

