

IV Armboards from Morrison Medical

Disposable IV Armboards

- Strong double wall corrugated core
- Comfortable foam padding
- Vinyl meets FDA requirements
- Dielectrically sealed
- Washable satin finish

Item #

1002-50	1 x 4", case of 50
1003-50	1½ x 3", case of 50
1004-50	2 x 4", case of 50
1005-50	2 x 6", case of 50
1005-100	2 x 6", case of 100
1010-50	2 x 9", case of 50
1010-100	2 x 9", case of 100
1008-50	3 x 6", case of 50
1009-50	3 x 9", case of 50
1009-100	3 x 9", case of 100
1012-50	3 x 12", case of 50
1012-100	3 x 12", case of 100
1016-50	3 x 15", case of 50
1016-100	3 x 15", case of 50
1019-50	3 x 18", case of 50
1019-100	3 x 18", case of 100



Reusable IV Armboards

- Strong wood core
- Comfortable foam padding
- Vinyl meets FDA requirements
- Dielectrically sealed
- Washable satin finish
- Premium quality, lightweight and affordable
- "L" shaped 4 x 5" footboards (padded on both sides) available to prevent infant foot drop

Item #

1032-50	1 x 4", case of 50
1020-50	1½ x 3", case of 50
1031-50	2 x 4", case of 50
1023-50	2 x 6", case of 50
1029-50	2 x 9", case of 50
1039-50	3 x 6", case of 50
1024-50	3 x 9", case of 50
1024-100	3 x 9", case of 100
1026-50	3 x 12", case of 50
1041-50	3 x 15", case of 50
1042-50	3 x 18", case of 50
1042-100	3 x 18", case of 100

For pricing and additional information, call 800-438-6677 or visit www.MorrisonMed.com.

