

## Economy Deluxe Limb Restraints & Quick Check Limb Holders

Morrison Medical Limb Restraints prevent combative patients from interfering with administration of medical care and harming themselves or medical personnel.

### Economy Deluxe Limb Restraints

- Constructed of soft breathable foam and strong 1" polypropylene strap.
- Velcro® closure holds restraint in place to allow quick application and tying of straps.
- Reusable, but economically priced to be considered disposable.
- Double 'D' ring style is a 11.5"x2.5" Cuff with single 60" strap.
- Tie Down style is a 11.5"x2.5" Cuff with two 1"x60" straps (adult model).
- Adult, pediatric and infant sizes available.
- Latex-Free.
- FDA Registered.

### Quick Check Limb Holders

- Constructed of soft breathable foam and strong 1" polypropylene strap.
- Velcro® closure holds restraint in place during application.
- Quick release buckle secures cuff to patient and allows for adjustment.
- Reusable, but economically priced to be considered disposable.
- Large 3" x 12" long cuff.
- Quick Check has a 60" single strap for tie down.
- Double Check has a 60" strap with a quick release buckle.
- Won't cinch like other restraints.
- Latex-Free.
- FDA Registered.

Item	1103	Adult Deluxe, tie down style (A), Pair
	1103-24	Adult Deluxe, tie down style (A), 24 Pair/Case
	1103-48	Adult Deluxe, tie down style (A), 48 Pair/Case
	1114	Adult Double D-ring style (B), Pair
	1114-48	Adult Double D-ring style (B), 48 Pair/Case
	1120	Infant Deluxe, tie down style, pair
	1120-48	Infant Deluxe, tie down style, 48 Pair/Case
	1124	Pediatric Deluxe, tie down style, pair
	1124-48	Pediatric Deluxe, tie down style, 48 Pair/Case
	1130	Quick Check limb restraint with buckle, Pair
	1132	Double Check limb restraint with buckle, Pair

Please call for pricing.



Deluxe, Tie Down Style (A)



Deluxe, D-Ring Style (B)



Quick Check, Item 1130



Double Check, Item 1132

For more information, call 800-438-6677,  
or visit [www.MorrisonMed.com](http://www.MorrisonMed.com).