

The Linebacker™ Tapeless Intravenous Tube Restraint

The Linebacker™ is easy to use, affordable, completely disposable and an effective way to secure a catheter tube to a body part. The product utilizes the simplicity of a comfortable strap secured by Velcro® that can be adjusted to custom fit patients. This is done in four easy steps:

Step 1: Intiate IV

Step 2: Peel back protective film on the Linebacker™

Step 3: Place IV tubing on adhesive and fold over

Step 4: Wrap strap around limb and secure



What makes The Linebacker™ the superior choice?

- Latex free
- Painless removal
- Quick and easy to apply
- Stronger, more secure hold than tape
- No potential of skin infection
- Can be used on hand, arm or foot
- Single unit can be used to secure tube as well as arm to IV armboard
- Single patient use, disposable
- Could be modified to secure other catheters such as Intra-aortic balloon pumps
- Works well on arterial lines as well as IVs

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*

