

Morrison Medical's Original Head Vise I



Morrison Medical's Original Head Vise I – the head immobilizer that started the Head Vise family of disposable and reusable head immobilizers. Made from heavy duty, impervious corrugated plastic, The Original Head Vise I is completely waterproof and can be written on for easy identification. It also is x-ray translucent, easy to disinfect for reuse, and stores flat to fit in compact spaces. This head immobilizer features a quick and sturdy application. The rugged exterior is lined with a comfortable foam pad for patient comfort. The Original Head Vise I fits virtually all patients with an adjustable feature to fit children and adults. It is also compatible with helmets. The Original Head Vise I is sold individually or in cases of 18. Straps and replacement foam inserts are also available.



*Folds flat for compact storage.
Assembles easily for quick use.*

Original Head Vise I *Please Call for Pricing*

Item #1450:
Original Head Vise I, each

Item #1450S:
Original Head Vise I w/ #1475 strap, each

Item #1450-18:
Original Head Vise I, case of 18

Item #1456-12:
Original Head Vise I replacement foam,
pack of 12

- White, heavy duty, impervious corrugated plastic.
- Completely waterproof with a writeable surface.
- Includes foam pad for patient comfort.
- X-ray translucent.
- Stores flat.
- Easily disinfected for reuse.
- Adjustable for children and adults.
- Compatible with helmets.
- Attaches to board quickly with EZ pull tabs.
- Works great with self adhesive disposable head/chin straps (#1475).

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*

