

Morrison Medical's Original Head Vise II

Morrison Medical's Original Head Vise II is made of sturdy, moisture resistant cardboard and features a wrap around foam pad for patient comfort. Features of The Original Head Vise II include: a writable surface for easy identification; x-ray translucent; stores flat to fit in compact spaces; and it is disposable. This head immobilizer features a quick and sturdy application. The Original Head Vise II fits virtually all patients with an adjustable feature to fit children and adults. The Original Head Vise II is sold individually or in cases of 10 and 50.



Original Head Vise II *Please Call for Pricing*

Item #1460:
Original Head Vise II, each

Item #1460S:
Original Head Vise II w/ #1475 strap, each

Item #1460-10:
Original Head Vise I, case of 10

Item #1460-50:
Original Head Vise I, case of 50

- Made of sturdy, moisture resistant cardboard.
- Includes comfortable foam head pad.
- X-ray translucent.
- Stores flat.
- One time use, disposable.
- Adjustable for children and adults.
- Attaches to board quickly with EZ pull tabs.
- Works great with self adhesive disposable head/chin straps (#1475).
- Writable surface.

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*