

Introducing The Complete Patient Transport System

Your confined space and bariatric solution.

Backboards

- Central concavity to center patient and decrease pressure and pain.
- Works as a pediatric board with special slotted area that uses standard straps to fit smaller patients. (Pediatric pins available)
- 100% x-ray translucent.
- 450, 600 and 1000 pound load capacity.
- Tapered for easy extrication in confined areas.
- Stores easily in emergency vehicles at 72" x 16" x 1¾".
- Custom monogramming available.
- Available in Red, Orange, Yellow, Blue, Green and Lime Green.


Load Loks

When using a backboard in combination with The Soft Stretcher, add extra security and patient support with Load Loks. Available in sets of 4, Load Loks stabilize the backboard to work cohesively with The Soft Stretcher. Velcro® closures allow for a custom fit for all board sizes, including pediatric.

The Soft Stretcher

Available in 3 sizes, Morrison Medical's Soft Stretchers are specially designed for moving patients of all sizes up to 2,000 pounds. Great for transporting patients out of awkward and confined areas or stairways, transferring patients, and moving several children or more than one average size adult at a time.

- 2,000+ pound load capacity
- Virtually impervious
- Latex Free
- Durable vinyl material with continuous reinforced nylon straps
- No-slip grip rubberized and pre-molded handles allow for safer transport and prevent transporter fatigue.
- Cleans quickly with normal cleaning protocol.
- Affordable – durable and long lasting means fewer necessary replacements.



Call for Pricing

#3010: Regular Soft Stretcher, 30" x 70"

#3020: Large Soft Stretcher, 40" x 71"

#3030: XXL Soft Stretcher, 60" x 80"

#3011: Regular Load Loks, Set of 4

#3021: Large Load Loks, Set of 4

#3031: XXL Load Loks, Set of 4

#3053: Backboard 450 Without Pins

#3054: Backboard 450 With Pins

#3055: Backboard 600 Without Pins

#3056: Backboard 600 With Pins

#3057: Backboard 1000 Without Pins

#3058: Backboard 1000 With Pins

*For more information, call 800-438-6677,
or visit www.MorrisonMed.com.*